

STUDIO



Shown in
Santa Barbara
(Fog)

Specifications	Standard	Petite
Weight Capacity	158.7kg	136kg
Seat Width	57.1cm	50.8cm
Seat Depth	54.6cm	45.7cm
Seat to Floor	46.9cm	45.7cm
Chair Width	96.5cm	91.4cm
Chair Depth	106.7cm	96.5cm
Chair Height	106.7cm	101.6cm
Operation	USB Hand Wand	USB Hand Wand

Modern design and incredible comfort, the Studio lift chair is sure to become a timeless classic. This space-saving dual motor Power Lift Recliner delivers a relaxation experience that's second to none.

Benefiting from an independent headrest, footrest extension, and power lumbar support, the Studio is available in three distinctive colour options.

Colour Options

Santa Barbara (Fog)

Santa Barbara (Latte)

Santa Barbara (Bark)

- Independent headrest
- Compact design
- Power lumbar support
- Footrest extension
- Choose from three colours
- USB-compatible hand control