

Fact Sheet



Bathroom



Shower Chair detachable Back

The Affinity Shower Chair with detachable backrest offers the user the option of stool or chair. Durable anodized aluminium frame with plastic blown seat and backrest. Backrest has multiple positioning options, can be removed to us as a stool. Rubber armrests provide comfort and safety.

- * Lightweight aluminium frame
- * Removable backrest
- * Drain holes to reduce water build up
- * Height adjustable
- * Large non-slip stoppers
- * Weight capacity 113kg



Specifications

Overall height.....	75-87	Seat height	41-54
Overall depth	43	Armrest height	15
Overall width.....	51	Unit weight.....	3.5
Seat depth.....	41-44	Weight capacity	113
Seat width	41		

Measurements in cm, weight in kg

Product Code: SCBAR

Your Local Agent

QLD 6-8 Imboon St, Deception Bay 4508
Ph 1300 368 085 qld@rehabandmobility.com.au
www.rehabandmobility.com.au