



Seating



Astral Perching Stool

The Astral Perching Stool has been developed for comfort and longevity, the contoured lines of both the seat and back are not only attractive to the eye, but also offer excellent ergonomic support with seamless polyurethane seat. The height adjustable frame has side mounted stress bars and a large footprint for safety and stability. Ideal for those who struggle to stand for long periods.

- * Ideal for perching when working in the kitchen or at a bench
- * Adjustable seat height
- * Padded backrest
- * Easy clean
- * Weight Capacity 159kg



Product Code: VG862



Product Code: VG863



Product Code: VG865

Specifications—VG865

Overall height..... 78-90
 Overall depth 41-65
 Overall width..... 46-65
 Seat depth.....30
 Seat width35
 Seat height 50-65

Backrest height27
 Width between arms48
 Armrest height17
 Unit weight.....5.2
 Weight capacity.....159

Measurements in cm, weight in kg

Your Local Agent

QLD 6-8 Imboon St, Deception Bay 4508
 Ph 1300 368 085 qld@rehabandmobility.com.au
www.rehabandmobility.com.au