

Fact Sheet

Scooter

Bandit

#1 Pull Apart



Australia's most comfortable new, compact, transportable mobility scooter. This new versatile yet lightweight Mobility Scooter has both front and rear suspension, easily disassembles for transport and storage into 5 lightweight pieces and has a powerful 300w motor.





FEATURES

- Easy dismantle for transport and storage
- 9" pneumatic tyres providing cushion from rough driving surfaces
- All around suspension for comfort
- Best power to weight ratio
- High visibility headlights for safety
- Swivel seat for easy mount and dismount
- Front to back and height adjustable seat
- Available colours: Blue, Red, Silver, White

Dismantled weights

- Seat 10.5kg
- Battery 13.5kg
- Transaxle 16.5kg
- Body 20.5kg

SPECIFICATIONS

Motor Power	300w/0.4hp
Max Speed	6.8km/h
Operation Range	Up to 25km
Max Climbing Grade	6-8° (depending on weight of user)
Turning Radius	1100mm
Ground Clearance	75mm
Max Loading Weight	136kg
Overall Dimensions	1080x505x870mm
Wheel Size– Front	9" (22cm) pneumatic
Weight with batteries	59.9kg
Battery	12V 20Ah x2 (lead-acid)
Controller	PG 45A
Brake System	Electromagnetic
Free Wheel Mode	Yes
Suspension	Front & Rear
Lights	Headlights
Charger Size	24V 2A

