



Fact Sheet

Walking Aids *Bariatric*

Bariatric Walking Frame

The Bariatric Walking Frame is a stable and helpful option for those needing additional support when walking. This frame is heavy duty, sturdy and height adjustable with a wide base. Foldable, making it easy for transportation and storage, bringing the sides in with two easy push buttons helps those with limited dexterity.



- * Sturdy heavy duty frame
- * Height adjustable
- * Easy push buttons to fold, transport and store
- * Dual front and side brace for extra support
- * Contoured foam handgrips for greater comfort
- * Can be customized with walker skis
- * Weight capacity 270kg

Walking Frame Accessories

VP179T Apron Bag

VP179N Net Bag

(details available www.rehabandmobilityshop.com.au)

Specifications

Overall height..... 80-105
Overall depth52
Overall width.....72
Folded width11

Unit weight.....5.7
Weight capacity.....270

Measurements in cm, weight in kg

Product Code: RM804

Your Local Agent

QLD 6-8 Imboon St, Deception Bay 4508
Ph 1300 368 085 qld@rehabandmobility.com.au
www.rehabandmobility.com.au