



Fact Sheet

Bathroom



Stackable Bath Step

The Plastic Bath Step features an extra-large platform with a textured, non-slip top surface and rubberised feet for increased safety. These handy bath steps can be securely stacked, raising the height in 51mm (2") increments, and allowing multiple units to be stacked and stored easily when not in use. The Plastic Bath Step is a versatile mobility aid that can also be used around the home, both inside and outside.



- * Versatile and lightweight
- * Secure, stackable step
- * Extra large platform
- * Non-slip textured surface
- * Rubber feet for safety
- * Weight capacity 120kg

Specifications

Overall height	10
Overall depth	42
Overall width	42
Unit weight	1
Weight capacity	120

Measurements in cm, weight in kg

Product Code: VR278D

Your Local Agent

QLD 6-8 Imboon St, Deception Bay 4508
Ph 1300 368 085 qld@rehabandmobility.com.au
www.rehabandmobility.com.au