

# Fact Sheet



## Bathroom



### Affinity Bath Transfer Bench

The Affinity Bath Transfer Bench is placed over the bath for those who have a shower over the bath situation, minimizing the risk of slips and falls. The wider seat allows the client to slide across from the outside of the bath. Four non-slip ferrules for secure placement.

- \* Aluminium frame
- \* Height adjustable
- \* Suction tips
- \* Easy clip together frame
- \* Weight capacity 180kg



#### Specifications

Overall height ..... 86-99  
 Overall depth ..... 75  
 Overall width ..... 53  
 Seat depth ..... 41  
 Seat width ..... 67

Tub clearance height . 41-53  
 Unit weight ..... 4.9  
 Weight capacity ..... 180

*Measurements in cm, weight in kg*

**Product Code: BA397HD**

**Your Local Agent**

**QLD** 6-8 Imboon St, Deception Bay 4508  
**Ph 1300 368 085** [qld@rehabandmobility.com.au](mailto:qld@rehabandmobility.com.au)  
[www.rehabandmobility.com.au](http://www.rehabandmobility.com.au)