

Fact Sheet



Bathroom *Bariatric*

Bath Transfer Bench

The durable aluminium Bath Transfer Bench is placed over the bath for those who have a shower over the bath situation. The wider seat allows the client to slide across from the outside of the bath. Two large suction feet for secure placement. Fully height adjustable will fit over most baths.

- * Durable aluminium frame
- * Detachable back rest, can be attached to either side
- * Heavy duty suction tips
- * Detachable arm rest
- * Weight capacity 200kg



- **Complies with Australian Standards
AS/NZS 3973:2009**

Specifications

Overall height.....	84-99	Tub clearance.....	41-55
Overall depth.....	55	Armrest height	30
Overall width.....	78	Unit weight.....	6
Seat depth	41	Weight capacity.....	200
Seat width	68		

Measurements in cm, weight in kg

Product Code: RM424

Your Local Agent

QLD 6-8 Imboon St, Deception Bay 4508
Ph 1300 368 085 qld@rehabandmobility.com.au
www.rehabandmobility.com.au