



Bathroom



Bewl Mobile Commode

The Bewl Shower Commode Chair boasts a rust-free, all alloy frame with a 25 year limited warranty, a seamless PU seat and back for maximum comfort and safety. The self propelled wheels on the Bewl Shower Commode require little or no maintenance. The Bewl Self Propelled Shower Commode Chair can easily be converted to an attendant propelled shower chair if the client's needs change.

- * Powder coated alloy frame
- * Self propelled
- * Contoured PU foam seat and backrest
- * 2 in 1 front cut away seat with solid infill
- * Quick release wheels
- * 2 locking 5" closed axis castors, resistant to hair & dirt
- * Foam padded flip back armrest
- * Adjustable, swing away, removable footplates
- * Pan with lid
- * Easy clean
- * Max toilet bowl height 435mm
- * Weight capacity 127kg

Also available as attendant propelled (VB503)

- * Max toilet bowl height 415mm



Specifications

Overall height	94	Seat height.....	52
Overall depth	108	Armrest height	18
Overall width	56	Unit weight	20
Seat depth	43	Weight capacity.....	127
Seat width.....	44		

Measurements in cm, weight in kg

Product Code: VB505

Your Local Agent

QLD 6-8 Imboon St, Deception Bay 4508
 Ph 1300 368 085 qld@rehabandmobility.com.au
 www.rehabandmobility.com.au