



Bathroom



Bewl Commode - Height Adjustable

The Bewl Adjustable Height Shower Chair is manufactured from high quality materials and boasts a rust free all alloy frame. As the name suggests the Bewl Shower Commode Chair is height adjustable allowing it to be set at the optimum height for the user. The seat and back support are manufactured from PU foam and are contoured for comfort and support. The height adjustable Bewl Shower Commode Chair also features a fold away footrest as standard offering a larger overall footprint for increased stability.



- * Powder coated alloy frame
- * Attendant propelled
- * Height adjustable
- * Contoured PU foam seat and backrest
- * 2 in 1 front cut away seat with solid infill
- * 4 locking 5" closed axis castors, resistant to hair & dirt
- * Foam padded flip back armrest
- * Adjustable, flip up footplate
- * Pan with lid
- * Easy clean
- * Weight capacity 127kg

Specifications

Overall height 99-104
 Overall depth 94
 Overall width 53
 Seat depth 44
 Seat width 44

Seat height 52-57
 Armrest height 23
 Unit weight 14
 Weight capacity 127

Measurements in cm, weight in kg

Product Code: VB503

Your Local Agent

QLD 6-8 Imboon St, Deception Bay 4508
 Ph 1300 368 085 qld@rehabandmobility.com.au
www.rehabandmobility.com.au