

Carbon Fibre Rollator

RC-HCT-9266A

.....



ROTHCARE

Features

.....

Lightweight carbon fibre frame
Large wheels for easy use indoors or outdoors
Handy storage bag & comfortable seat

Tech Specs

.....

Width inside handgrips	47cm
Overall width	64cm
Overall depth	69cm
Seat Area	46x25cm
Seat Height	54cm
Adjustable handle height	79.5-92cm
Front wheel	10 inch (PVC)
Rear wheel	8 inch (PVC)
Max. user weight	136kg
Weight	5.9kg



Your Local Agent

QLD 6-8 Imboon St, Deception Bay 4508
Ph 1300 368 085 qld@rehabandmobility.com.au
www.rehabandmobility.com.au