



# Fact Sheet

## Incontinence



### Chair Pad



**Product Code: VM842** (Brown)  
**VM842A** (Blue)

- \* Hygienic protection for sofas, chairs and beds
- \* Soft, breathable top layer with a non-slip, waterproof backing
- \* Special stay-dry fabric for longer-lasting dryness and comfort
- \* A more eco-friendly option than disposable pads
- \* Machine-washable for total convenience
- \* Dimensions 590 x 470mm
- \* Absorption 800ml

### Bed Pad with Tuckins

- \* A more eco-friendly option than disposable bed pads
- \* Attached 'tuck-ins' fold around the mattress for a secure fit
- \* Soft breathable top layer with a non-slip, waterproof backing
- \* Special stay-dry fabric for longer-lasting dryness and comfort
- \* Machine-washable for total convenience
- \* Dimensions 850x900x25
- \* Absorption 1000ml



**Product Code: VM842B**

#### Your Local Agent

**QLD** 6-8 Imboon St, Deception Bay 4508  
Ph 1300 368 085 [wholesale@rehabandmobility.com.au](mailto:wholesale@rehabandmobility.com.au)  
**NSW** Unit 2/144 Lake Entrance Road, Oak Flats 2529  
Ph 1300 002 805 [nsw@rehabandmobility.com.au](mailto:nsw@rehabandmobility.com.au)  
[www.rehabandmobility.com.au](http://www.rehabandmobility.com.au)