



Walking Aids

Compact Indoor Rollator

This lightweight sturdy compact rollator with a height adjustable handle is perfect for indoor and around the home use. With a delta handle and brake it suits those with limited upper limb mobility, stroke survivors or limb amputees. The removable tray and storage bag makes it easy to transport meals comfortably, great for carrying utensils or personal belongings around the house. This rollator is easy to fold, ideal for storage and transporting.

- * Lightweight sturdy frame
- * Delta handle and brake
- * Perfect for indoor use around the home
- * Handy tray for carrying meals
- * Removable storage bag
- * 8" Castors
- * Easy folding for storage and transport
- * Adjustable handle height



Specifications

Overall height85-96
 Overall depth 75
 Overall width 58
 Folded width 42

Tray size32x23x2
 Unit weight 6.7
 Weight capacity 136
Measurements in cm, weight in kg

Product Code: RM926G

Your Local Agent

QLD 6-8 Imboon St, Deception Bay 4508
 Ph 1300 368 085 qld@rehabandmobility.com.au
www.rehabandmobility.com.au