



# Fact Sheet

## Pressure Care



## Cool Comfort Cushion

- \* Ideal for use in the home, car or as a wheelchair cushion
- \* Ergonomically designed memory foam contoured to improve posture & comfort
- \* Absorbs up to 84% of body pressure
- \* Unique vent hole design to aid cooling
- \* Cover offers

Cooling fabric top reduces body temp by 2°

Breathable mesh sides

Anti-slip base

- \* Available in 3 sizes
- \* Weight Capacity 220kg



### Cool Comfort Pressure Absorbing Cushion

The Affinity Cool Comfort Cushion is designed to absorb body pressure and reduce back stress. The Ergonomic design is comfortable, and encourages good posture.

The unique combination of cooling fabric cover, vent design and breathable mesh side fabric makes the Affinity Cool Comfort Cushion perfect for long term, everyday, use.

Absorbs up to 84%  
of Body Pressure & Reduces  
Body Temperature by 2°

Unique Technology  
Ergonomic Design  
SGS, ROHS, & CE Certified  
For your Peace of Mind



Suitable for Seat & Wheelchair

**Washing Instructions**  
1) Remove cover from cushion and wash separately  
2) Use mild detergent  
3) Dry on a flat surface

Top: Delicate Cooling Fabric

Sides: Breathable Mesh

Base: Anti-slip Backing



### Contoured & Vented Design





Unique Vent-Hole  
Design to Aid Cooling



Ergonomically Designed to  
Improve Posture & Comfort

### Pressure Mapping Data



## Specifications

Product Code: CSH-CCL 44x42x8cm

CSH-CCM 42x40x8cm

CSH-CCC 50x46x6cm

[www.affinityproducts.com.au](http://www.affinityproducts.com.au)



**Your Local Agent**

QLD 6-8 Imboon St, Deception Bay 4508  
Ph 1300 368 085 [wholesale@rehabandmobility.com.au](mailto:wholesale@rehabandmobility.com.au)  
NSW Unit 2/144 Lake Entrance Road, Oak Flats 2529  
Ph 1300 002 805 [nsw@rehabandmobility.com.au](mailto:nsw@rehabandmobility.com.au)  
[www.rehabandmobility.com.au](http://www.rehabandmobility.com.au)