



Fact Sheet

Wheelchair



Deluxe Self Propelled

The Deluxe Lightweight Self Propelled Aluminium Wheelchair from Aidapt is our latest design in lightweight, compact wheelchairs. High grade aluminium with a double cross bar for strength, the lightweight wheelchair also offers ease of movement for the user and quick-folding for trouble-free transportation. Ideal for the active, independent user.

- * Lightweight Self Propelled
- * Care assist hand brake
- * Folding backrest (for storage or transport)
- * Quick release wheels (for transport)
- * Flip up armrest for side transfer
- * 24" rear wheel, 8" front castor
- * Swing away and removable leg rests
- * Solid tyres (no maintenance required)
- * Can be fitted with elevating leg rests
- * Weight capacity 115kg



Specifications

Overall height 89
 Overall depth 108
 Overall width 63
 Seat depth 40
 Seat width 46

Seat height 51
 Backrest height 40
 Armrest height 23
 Unit weight 14
 Weight capacity 115

Measurements in cm, weight

Product Code: VA165B Blue
Product Code: VA165R Red
Product Code: VA165S Silver
Product Code: VA165O Orange

Your Local Agent

QLD 6-8 Imboon St, Deception Bay 4508
Ph 1300 368 085 qld@rehabandmobility.com.au
www.rehabandmobility.com.au