

## Walking Aids

### Duo Transporter Rollator

The Duo Transporter is a high quality mobility aid having a dual purpose, converting from a rollator to a wheelchair in seconds. This allows the user to walk with support of the rollator and to sit down and rest when needed. The adjustable, removable foot rests allow it to be used as a wheelchair with a flip able backrest to enable the user to be pushed. Suitable for indoor and outdoor use it is provided with a motion and park brake, and height adjustable handle bars. The strong aluminium frame is easily foldable and great for frequent travellers.



Product Code: RM199GD

- \* Aluminium frame
- \* Double cross brace frame
- \* Great for frequent travellers
- \* Height adjustable handle bars
- \* Motion and park brake
- \* Flip able padded backrest
- \* Easy to fold and transport
- \* Removable bag
- \* Max weight user: 110kg
- \* Colour: Gold

#### Specifications

Overall Height	Overall Length	Overall Width	Seat Depth	Seat Width	Seat Height	Folded Width	Unit Weight	Weight Capacity
85-97	Rollator:74 Wheelchair:101	63	30	48	55.5	28	9.7	110

Measurements in cm, Weight in kg

**Your Local Agent**

**QLD** 6-8 Imboon St, Deception Bay 4508  
 Ph 1300 368 085 [wholesale@rehabandmobility.com.au](mailto:wholesale@rehabandmobility.com.au)  
**NSW** Unit 2/144 Lake Entrance Road, Oak Flats 2529  
 Ph 1300 002 805 [nsw@rehabandmobility.com.au](mailto:nsw@rehabandmobility.com.au)  
[www.rehabandmobility.com.au](http://www.rehabandmobility.com.au)