

Fact Sheet



Bedroom



Essex Bedside Commode

The Essex Commode Chair is height adjustable, practical commode chair suitable for use as a bedroom commode or bathroom commode. The Essex Commode Chair has a 5 litre bucket with lid and a padded over seat is also supplied for discretion not in use. The upholstery of the commode is well padded and is vinyl welded to reduce the possibility of bacteria build up. Available fully assembled or flat packed.



- * Padded back
- * Detachable padded over seat
- * Height adjustable
- * 5 litre pail with lid supplied
- * Hinged thermoset plastic toilet seat
- * Flat pack for transportation and storage
- * Weight Capacity 158kg

Specifications

Overall height 82-103
 Overall depth 46
 Overall width 54
 Seat depth 42
 Seat width 42

Seat height 49-70
 Armrest height 17
 Unit weight 6.3
 Weight capacity 158

Measurements in cm, weight in kg

Product Code: VR161G

Your Local Agent

QLD 6-8 Imboon St, Deception Bay 4508
Ph 1300 368 085 qld@rehabandmobility.com.au
www.rehabandmobility.com.au