

Fact Sheet

Scooter

Everest

The #1 Choice



The Everest cruises up difficult inclines with ease and navigates over rough terrain with clearance to spare. Dominate the streets with a class leading turning circle and a powerful high torque motor. This robust go anywhere scooter has an unequalled array of safety and comfort features. With an amazing load carrying capacity the Everest is streets ahead of the competition and commands respect with its proud looks and modern finish.





FEATURES

- Industry leading high torque motor and gearbox
- Safety turning sensor
- USB charging point
- 203kg weight capacity
- Super comfort pillow to seat with high back
- Rider comfort seat post suspension
- High ground clearance
- Safety park brake
- Digital dash display
- Industry leading turning circle
- Independent front suspension
- Front and rear adjustable suspension
- Available colours: Orange, Silver



SPECIFICATIONS

Motor Power	1000w/1.3hp
Max Speed	10km/h
Operation Range	Up to 49km
Max Climbing Grade	12° (depending on weight of user)
Turning Radius	1600mm
Ground Clearance	195mm
Max Loading Weight	203kg
Overall Dimensions	1450x710x1342mm
Wheel Size– Front	13" (33cm) pneumatic
Weight with batteries	146kg
Battery	12V 75Ah x2 (lead-acid)
Controller	PG S-140A
Brake System	Electromagnetic 8Nm
Free Wheel Mode	Yes
Suspension	Front & Rear
Rear Mirrors	Left & Right
Charger Size	24V 8A

