



## Walking Aids

### Forearm Walker

The aluminium lightweight forearm rollator allows the user to support their weight on the arms cuffs to reduce strains on the wrists and hands. The forearm cuffs are contoured and height adjustable. Featuring 10" front and 8" rear PVC wheels, and includes a detachable shopping bag. Easily foldable and compact for storage.



- \* Aluminium lightweight frame
- \* Height adjustable handles
- \* Easily foldable and compact
- \* Includes fabric seat and detachable shopping bag
- \* 10" front and 8" rear PVC wheels
- \* Forearm gutter height 94-109cm
- \* Weight capacity 136kg

#### Specifications

Overall height.....	106-121	Seat height .....	55
Overall depth.....	67	Armrest height .....	94-109
Overall width.....	62	Folded width .....	29
Seat depth .....	23	Unit weight.....	12
Seat width .....	44	Weight capacity.....	136

*Measurements in cm, weight in kg*

**Product Code: VP187F**

**Your Local Agent**

**QLD** 6-8 Imboon St, Deception Bay 4508  
**Ph 1300 368 085** [qld@rehabandmobility.com.au](mailto:qld@rehabandmobility.com.au)  
[www.rehabandmobility.com.au](http://www.rehabandmobility.com.au)