

# Fact Sheet



## Wheelchair

### G3 Wheelchair

The G3 Wheelchair is equipped with a number of special features along with a host of accessory choices, makes this a very desirable chair for everyday use. Care assist brake offering reassurance to the carer when on ramps or inclines. 12" rear castor for a more comfortable ride.



- \* Sturdy aluminium frame
- \* Angle and height adjustable seat
- \* Folding backrest
- \* Care assist brake
- \* Swing away removable leg rests
- \* Padded nylon upholstery
- \* 12" rear castor
- \* 7" front castor
- \* Anti-tip wheels
- \* Weight capacity 150kg



Optional Stump  
Support



Optional Elevating  
Leg Rest

### Specifications

Overall height .....	75-90	Seat height.....	46-50
Overall depth .....	107	Backrest height .....	38
Overall width .....	67	Armrest height.....	22
Seat depth.....	41	Unit weight .....	16
Seat width.....	46	Weight capacity .....	150

*Measurements in cm, weight*

**Product Code: RM138-46T**

**Your Local Agent**

**QLD** 6-8 Imboon St, Deception Bay 4508  
**Ph 1300 368 085** [qld@rehabandmobility.com.au](mailto:qld@rehabandmobility.com.au)  
[www.rehabandmobility.com.au](http://www.rehabandmobility.com.au)