



Wheelchair

G3 Wheelchair

The G3 Wheelchair is equipped with a number of special features along with a host of accessory choices, makes this a very desirable chair for everyday use. The self propel model is available with a standard brake or a drum brake along with care assist handbrake.



- * Sturdy aluminium frame
- * Angle and height adjustable seat
- * Folding backrest
- * Care assist brake
- * Quick release 24" mag rear wheels with PU Tyres
- * Swing away removable leg rests
- * Padded nylon upholstery
- * 7" front castor
- * Optional drum brake (S/P only)
- * Anti-tip wheels
- * Weight capacity 150kg

Optional Stump
Support



Optional Elevating
Leg Rest



Specifications 46cm/51cm

Overall height 74-90
Overall depth 107
Overall width 67/71
Seat depth 41
Seat width 46/51

Seat height 46-50
Backrest height 38
Armrest height 22
Unit weight 17
Weight capacity 150

Measurements in cm, weight

Product Code: RM138-46 46cm Seat
Product Code: RM138-46D Drum Brake
Product Code: RM138-51 51cm Seat
Product Code: RM138-51D Drum Brake

Your Local Agent

QLD 6-8 Imboon St, Deception Bay 4508
Ph 1300 368 085 qld@rehabandmobility.com.au
www.rehabandmobility.com.au