



Fact Sheet

Wheelchair

G6 Excel Bariatric Wheelchair

The G6 Excel Bariatric Wheelchair is multi-adjustable with individual tailored options allowing this chair to be modified to the individual user. A number of important features along with a host of accessory choices, makes this a very desirable chair for everyday use with the comfort of a padded seat. Care assist brake gives a little extra assurance for the career. This model is available with a standard brake or with a drum brake providing extra comfort for the passenger.



- * Strong steel frame
- * Angle and height adjustable seat
- * Swing away removable leg rests
- * Flip back and removable armrests
- * Quick release 24" rear wheels
- * 7" front caster
- * Attendant assist hand brake
- * Optional drum brake
- * Reclining backrest
- * Weight capacity 250kg



Optional Stump Support



Optional Elevating Leg Rest



Optional 1 Piece Footplate

Specifications 56/61

Overall height103-116
 Overall depth 102
 Overall width 78/82
 Seat depth..... 46
 Seat width 56/61

Seat height..... 53
 Backrest height 45
 Armrest height..... 23
 Unit weight 26/27
 Weight capacity 250

Measurements in cm, weight

Product Code: RM160-56 56cm Seat
Product Code: RM160-56D Drum Brake
Product Code: RM160-61 61cm Seat
Product Code: RM160-61D Drum Brake

Your Local Agent

QLD 6-8 Imboon St, Deception Bay 4508
 Ph 1300 368 085 qld@rehabandmobility.com.au
www.rehabandmobility.com.au