



Bathroom



Grab Rails

Grab rails are installed to provide greater independence, providing stability when using bathroom facilities and for general use around the home.

Solo Contract Hinged Arm Support

Product Code: VY423



- * Ideal for use in the bathroom where safety is a concern
- * Simple fold up / fold down mechanism with arm support leg
- * Epoxy coated ensuring durability
- * 750mm long bar

Specifications

Depth to Wall (mm).....	760
Width (mm)	100
Leg Height (mm).....	710-870
Leg from wall (mm)	555
Weight Capacity (kg)	115

Solo Easigrip

- * Durable epoxy coated steel
- * Rounded ends for safety



Product Code: 300mm (12") VY393



450mm (18") VY394

President Ribbed Grab Bar

- * Durable UPVC finished in white
- * Ribbed grip for improved safety
- * Reinforced fixing points with cover caps
- * High quality AVS flanges
- * 32mm tube
- * Depth from wall 60mm



Product Code: 300mm (12") VY448



Product Code: 300mm (12") VY451
450mm (18") VY452



Product Code: 600mm (24") VY427



Product Code: 450mm (18") VY426



Product Code: 300mm (12") VY425

Your Local Agent

QLD 6-8 Imboon St, Deception Bay 4508
Ph 1300 368 085 wholesale@rehabandmobility.com.au
NSW Unit 2/144 Lake Entrance Road, Oak Flats 2529
Ph 1300 002 805 nsw@rehabandmobility.com.au
www.rehabandmobility.com.au