



Bathroom



Hythe Mobile Shower Chair

The Hythe Mobile Shower Chair incorporates the latest in seat design for improved safety and comfort. Height adjustable with 4 locking castors and features a pram handle as standard. The seat drainage outlets fully conform to current MHRA safety guidelines. Available with a variety of options for height adjustment and brakes.

- * Tubular steel frame with tough polymer coating
- * Pram handle for carer
- * Moulded polypropylene seat
- * 4 large lock in 5" castors
- * Easy to clean and maintain
- * Mobile shower chair ideal for assisting users in and out of the shower area
- * Can be used in wet environments and level access showers
- * Weight capacity 127kg



Specifications

Overall height 83-92
 Overall depth 63
 Overall width 49
 Seat depth 40
 Seat width 44

Seat height 51-60
 Armrest height 15
 Unit weight 7
 Weight capacity 127

Measurements in cm, weight in kg

Product Code: VB650

Your Local Agent

QLD 6-8 Imboon St, Deception Bay 4508
 Ph 1300 368 085 qld@rehabandmobility.com.au
 www.rehabandmobility.com.au