



Walking Aids

Aidapt Knee Walker

The Innovative Knee Walker is a great option to help with the recovery from a foot injury, break or sprain, allowing the user to evenly distribute their body weight on the device. With the supplied cushioned knee pad to support injured leg, ankle or foot providing maximum comfort. It has a sturdy dual frame which adds stability and durability helping with easy steering and stopping.

- * Sturdy dual frame
- * Compact for easy transportation
- * Height adjustable handle bar and knee pad
- * Detachable basket
- * Fits adults, young adults and elderly.
- * Weight capacity 113kg



Additional Sheepskin Knee Pad Cover available

Product Code: SS6000



Specifications

Overall height86-106
 Overall depth 80
 Overall width 43
 Seat depth..... 8
 Seat width..... 70

Seat height.....48-58
 Unit weight 10
 Weight capacity 113

Measurements in cm, weight in kg

Product Code: VP186

Your Local Agent

QLD 6-8 Imboon St, Deception Bay 4508
 Ph 1300 368 085 qld@rehabandmobility.com.au
www.rehabandmobility.com.au