



Seating Bariatric

Low Back Chair

The Low Back Orthopaedic chair is perfect for post therapy recovery. With both height adjustable and seat depth adjustment, this chair is suited to most users. Works perfectly as a dining chair with rubber stoppers to protect floors. Easy clean vinyl makes it ideal for both home or facility use. Available in two seat widths and weight capacities.



- * Ideal for post therapy comfort
- * Strong sturdy construction
- * Powder coated steel frame,
- * Padded seat & armrests
- * Height adjustable
- * Dismantle for easy transportation
- * Easy clean Vinyl
- * Available in Grey
- * Weight Capacity 205kg
Bariatric 250kg



Specifications

Overall height..... 42-57
 Overall depth59
 Overall width.....37
 Unit weight.....4.3

Product Code: RM531LG

Specifications—Standard/Bariatric

Overall height.....	89-105	Seat height	48-65
Overall depth	49/49	Armrest height	15
Overall width.....	55/69	Unit weight.....	12/13
Seat depth.....	43-46	Weight capacity.....	205/250
Seat width	46/60		

Measurements in cm, weight in kg

Product Code: RM535 Standard
Product Code: RM535B Bariatric

Your Local Agent

QLD 6-8 Imboon St, Deception Bay 4508
Ph 1300 368 085 qld@rehabandmobility.com.au
www.rehabandmobility.com.au