



Fact Sheet

Bedroom

Memory Foam Mattress

This 3 core foam mattress is made up of 3 different high quality foam types. The base layer is a 35D grade foam, support layer is a 30D grade foam and the top layer being a 60D grade heavy duty memory foam. This mattress maximizes immersion and envelopment of body to achieve maximum contact. Providing comfort and passive pressure relief protecting against ulcers and pressure sores. The foams are viscoelastic, free from latex, temperature sensitive and moulds to users shape.



- * 3 layers of foam (latex free)
- * 4 way stretch PU cover
- * Waterfall zip on 3 sides for easy removal
- * Machine washable up to 95degrees
- * Fire retardant foam and cover
- * Single: 1980x880x140mm
- * King single available on request: 1980x1050x150mm
- * Recommended patient weight 30-150kg

Product Code: M02-S

Your Local Agent

QLD 6-8 Imboon St, Deception Bay 4508
Ph 1300 368 085 qld@rehabandmobility.com.au
www.rehabandmobility.com.au