

Fact Sheet



Bathroom



Medina Bath Board

The Medina Bath Board is a high quality, reinforced plastic bath board.

- * Modern design
- * Slatted surface provides adequate drainage
- * Wipe clean surface
- * Simple to install with self locking cleats
- * Both available with or without handle
- * Weight Capacity 159kg

Specifications

Overall depth	4 slat 250
Overall depth	6 slat 350
Seat height	140
Seat width	685 / 710
Width adjustment	515-675 / 510-670
Net weight.....	3.2kg

Measurements in mm, weight in kg



Medina Bath Seat—Bariatric

The Medina Bath Seat provides extra height enabling easy access in and out of the bath.

- * Constructed from high quality, fully reinforced slatted plastic
- * Available in 2 heights (6" & 8")
- * Suitable for majority of bath's
- * Slatted for easy drainage
- * Modern easy to clean surface
- * Secured to bath with 4 large suction caps
- * Easily removed by pulling on the integral tabs
- * Finished in soft grey
- * Weight Capacity 190kg

Specifications

Dimensions (6") D290xW450xH150

Dimensions (8") D290xW450xH200

Net weight 2.4kg



Product Code: VR127A (6")

VR128A (8")

Your Local Agent

QLD 6-8 Imboon St, Deception Bay 4508
Ph 1300 368 085 qld@rehabandmobility.com.au
www.rehabandmobility.com.au