



Bathroom



Milton Aluminium Shower Chair

The Milton Stackable Shower Chair is a portable shower chair providing secure seat for those unable to stand in the shower. Incorporating the latest design in moulded shower seat for improved safety, drainage and comfort. Fixed armrests allow for safe transfer in and out of the shower.



- * Lightweight rust proof aluminium frame
- * Durable powder coat finish
- * Height adjustable
- * Moulded plastic seat with drain holes
- * Non-slip rubber tips
- * Weight capacity 150kg

• **Australian Standards**
AS/NZS 3973:2009

Specifications

Overall height..... 79-91
Overall depth44
Overall width.....57
Seat depth.....41
Seat width44

Seat height42-55
Armrest height 20
Unit weight.....3.7
Weight capacity 150

Measurements in cm, weight in kg

Product Code: VB499AL

Your Local Agent

QLD 6-8 Imboon St, Deception Bay 4508
Ph 1300 368 085 qld@rehabandmobility.com.au
www.rehabandmobility.com.au