

Scooter

Mustang

Luxury, Power, Comfort



With an array of safety features, all round double sprung suspension, seat post suspension, a luxury seat with adjustable lumbar support and multiple driver adjustment options makes this ride very comfortable. The Mustang has the most powerful motor and battery combination on the Top Gun range





FEATURES

- High performance 1650W 4 pole motor
- Safety turning sensor
- USB charging port
- 180kg weight capacity
- Luxury high back seat with lumbar support
- Rider comfort seat post suspension
- Dual tiller operation
- Safety park brake
- Digital dash display
- Front and rear double sprung suspension
- Available colours: Champagne, Red



SPECIFICATIONS

Motor Power	1650w/2.2hp
Max Speed	10km/h
Operation Range	45-60km
Max Climbing Grade	12° (depending on weight of user)
Turning Radius	2000mm
Ground Clearance	145mm at mid point
Max Loading Weight	180kg
Overall Dimensions	1461x690x1317mm
Wheel Size– Front	13" (33cm) pneumatic
Wheel Size—Rear	14" (35cm) pneumatic
Weight with batteries	150kg
Battery	12V 100Ah x2 (lead-acid)
Controller	PG S-200A
Brake System	Electromagnetic 10Nm
Free Wheel Mode	Yes
Suspension	Front & Rear
Rear Mirrors	Left & Right
Charger Size	24V 8A



Live a life with freedom

