



## Wheelchair



### Neptune 16" Wheelchair

The Neptune Self-Propel Wheelchair is a great option for everyday use with the large pneumatic self propel wheels giving a smooth ride. Featuring both a folding frame and backrest making it easy to store and transport into vehicles. With swing away, removeable leg rests and flip up arm rests this chair is great for the user to get in and out with ease, or when transferring sideways to and from the chair.

- \* Lightweight self propelled
- \* Care assist hand brake
- \* Folding backrest (for storage or transport)
- \* Quick release wheels (for transport)
- \* Flip up armrest for side transfer
- \* 24" rear wheel, 8" front castor
- \* Swing away and removable leg rests
- \* Weight capacity 100kg



#### Specifications

Overall height .....	92	Seat height.....	47
Overall depth .....	100	Backrest height .....	42
Overall width .....	57	Armrest height.....	20
Seat depth.....	40	Unit weight .....	13
Seat width.....	41	Weight capacity .....	100

*Measurements in cm, weight*

**Product Code: 9346376060721**

**Your Local Agent**

QLD 6-8 Imboon St, Deception Bay 4508  
 Ph 1300 368 085 qld@rehabandmobility.com.au  
 www.rehabandmobility.com.au