



## Bathroom

### Combo Shower Stool/Toileting Aid with Insert

The Combo Shower Stool/Toileting Aid offers an open gap in front and rear seat allowing access to body areas usually hard to reach. Can be used as a combination of a commode or provide a secure safe option for those who require a seat in the shower. Suitable for home use, hospitals and nursing homes. Black seat insert included.

- \* Lightweight rust proof aluminium frame
- \* Durable powder coat finish
- \* Height adjustable extension legs
- \* Suction rubber tips to reduce slip
- \* Soft polyurethane padded seat
- \* Includes bucket with lid
- \* Very user friendly for both the carer and user
- \* Black seat insert included
- \* Weight capacity 160kg



#### Specifications

Overall height ..... 65-81  
 Overall depth ..... 61  
 Overall width ..... 67  
 Seat depth ..... 45  
 Seat width ..... 48

Seat height ..... 42-57  
 Armrest height ..... 21  
 Unit weight ..... 6.3  
 Weight capacity ..... 160

*Measurements in cm, weight in kg*

**Product Code: RM370-48**

**Your Local Agent**

**QLD** 6-8 Imboon St, Deception Bay 4508  
 Ph 1300 368 085 [qld@rehabandmobility.com.au](mailto:qld@rehabandmobility.com.au)  
[www.rehabandmobility.com.au](http://www.rehabandmobility.com.au)