



Bathroom



Pembury Corner Stool

The Pembury height adjustable corner shower stool, an excellent option when space is at a premium. Featuring fine leg adjustment maximizing stability and safety when using the shower. A clip on/clip off seat with drainage holes allows for easy cleaning. Constructed from lightweight aluminium makes it easy to manage.



- * Ideal where space is at a premium
- * Powder coated
- * Height Adjustable
- * Non-slip rubber tips
- * Weight capacity 159kg

Specifications

Overall height 47-62
 Overall depth 38
 Overall width 52
 Seat depth 26
 Seat width 53

Seat height 47-62
 Unit weight 1.9
 Weight capacity 159

Measurements in cm, weight in kg

Product Code: VB636

Your Local Agent

QLD 6-8 Imboon St, Deception Bay 4508
 Ph 1300 368 085 qld@rehabandmobility.com.au
www.rehabandmobility.com.au