



## Wheelchair



## Ramps

### Lightweight Suitcase

- \* Suits majority of scooters and wheelchairs
- \* Provides access to vehicles and buildings
- \* Safety all-weather tread
- \* Easy to use and transport

#### Specifications

Folded width ..... 370mm  
 Channel width..... 700mm  
 Weight capacity 272kg (WCR180).....300kg



#### Product Code:

**VA142C** 910mm(3ft) net weight 9kg  
**VA142K** 1520mm(5ft)net weight 14kg  
**WCR180**1830mm(6ft) net weight 19kg

### Telescopic Channel

- \* Gritted surface for increased grip
- \* Lightweight

#### Specifications

Length collapsed (VA147F) ..... 744mm  
 (VA147S)...1201mm  
 Channel width ..... 189mm  
 Weight capacity ..... 180kg



#### Product Code:

**VA147F** 1220mm(4ft) net weight 6kg  
**VA147S** 2130mm(7ft) net weight 10kg

### Ramp Roll Up

Product Code: (3ft) VA141



- \* Easy access to houses and cars
- \* Portable & Compact
- \* Supplied with handy carry bag (20x20x90cm)

#### Specifications

Length ..... 760mm  
 Width .....200mm  
 Height .....900mm  
 Weight capacity ....300kg  
 Net weight..... 9kg

### Folding Suitcase

- \* Non slip surface
- \* Tri-fold for easy storage and transport



#### Specifications

Folded width ..... 360mm  
 Folded length .....630 &915mm  
 Channel width ..... 740mm  
 Weight capacity ..... 272kg

#### Product Code:

**VA143D** 1220mm(4ft) net weight 11kg  
**VA143S** 1830mm(6ft) net weight 20kg