



Fact Sheet

Walking Aids



Capero Lightweight x-fold Rollator

With the Capero aluminium lightweight rollator you will enjoy your new found freedom. Weighing just a little over 6kg makes the Capero a very affordable option in lightweight rollators. When in need of a break you have a comfortable nylon fabric seat and supportive backrest. Net bag and walking stick/umbrella holder complete this chic rollator.

- * Very Lightweight Aluminium frame
- * Easy side fold
- * Height adjustable handles
- * Canvas seat with padded backrest
- * 8" front castor 7" rear castor
- * Walking stick/umbrella holder
- * Detachable shopping bag
- * Net weight (no bag) 6kg
- * Max weight user: 136kg



Specifications

Overall height86-97
 Overall depth 76
 Overall width 64
 Seat depth..... 25
 Seat width..... 45

Seat height..... 54
 Folded width 25
 Unit weight 6.2
 Weight capacity 136

Measurements in cm, weight

Product Code: CAPEROMBL Blue
Product Code: CAPEROPLA Platinum

Your Local Agent

QLD 6-8 Imboon St, Deception Bay 4508
Ph 1300 368 085 qld@rehabandmobility.com.au
www.rehabandmobility.com.au