



# Fact Sheet

## Walking Aids



### Rollator Flexi Height Adjustable

The Affinity Flexi is light and nimble, with a unique seat height adjustment feature. Offering 20 height combinations making the Flexi the perfect custom fit for most users.

- \* Lightweight aluminium frame
- \* Metallic paint finish
- \* Height adjustable seat
- \* Ergonomic handles for comfort
- \* Fabric padded seat, fabric bag
- \* Walking stick holder
- \* Max weight user: 120kg



**Product Code:** 4W-6SHLG 6inch Castor  
Lime Green    4W-8SHALG 8inch Castor



**Product code:** 4W-6SHBO 6inch Castor  
Burnt Orange    4W-8SHAIO 8inch Castor



**Product Code:** 4W-6SHIB 6inch Castor  
Ice Blue    4W-8SHAIB 8inch Castor



**Product Code:** 4W-6SHIC 6inch Castor  
Charcoal    4W-8SHAIC 8inch Castor

#### Specifications 6"/8"

Overall height . 78-95/82-99  
Overall depth ..... 61/64  
Overall width ..... 61  
Seat depth..... 25  
Seat width ..... 38

Seat height..... 49-57/53-61  
Folded width ..... 23  
Unit weight ..... 6.5  
Weight capacity ..... 120

*Measurements in cm, weight*

**Your Local Agent**

**QLD** 6-8 Imboon St, Deception Bay 4508  
Ph 1300 368 085    qld@rehabandmobility.com.au  
[www.rehabandmobility.com.au](http://www.rehabandmobility.com.au)