



Fact Sheet

Bathroom

Rotating Seat Shower Chair

This Shower Chair features a rotating seat which rotates in both directions 360 degrees with lock positions at 90 degree intervals making it highly recommended for being user friendly to both the user and the carer. Constructed from a durable oval shaped aluminium frame with a detachable back and arm rests, height adjustable extension legs and suction rubber tips to resist slips.

- * Oval shaped aluminium frame
- * 360 degree rotating seat
- * Drainage holes
- * Removable backrest
- * Removable armrests
- * Height adjustable extension legs
- * Non-slip stoppers
- * Weight capacity 115kg



Specifications

Overall height..... 85-97
 Overall depth55
 Overall width.....58
 Seat depth.....33
 Seat width51

Seat height50-62
 Armrest height34
 Unit weight.....6
 Weight capacity115

Measurements in cm, weight in kg

Product Code: RM410

Your Local Agent

QLD 6-8 Imboon St, Deception Bay 4508
Ph 1300 368 085 qld@rehabandmobility.com.au
www.rehabandmobility.com.au