



Shifty Transfer Assist Trolley

The Shifty Transfer Assist Trolley with Leg Spread is designed to assist in transferring a person from one seated position to another. It is suitable for use in homes or healthcare facilities and features a small narrow frame that allows easy access through doors and corridors. It is smoothly powder coated with padded foam sections, both designed to be easily cleaned at high temperatures. The user is required to be able to weight bear through their legs and arms.



- * Suitable for in homes and healthcare facilities
- * Easy clean
- * Pedal operated leg spread
- * Swing away seat support
- * Support belt included
- * 2 locking castors
- * Weight capacity 200kg

Australian Standards AS/NZS ISO 10535-2011

Specifications

Overall height.....	107
Overall depth	80
Overall width.....	61.5-72.9
Unit weight.....	24
Weight capacity.....	200

Measurements in cm, weight in kg

Product Code: 9346376050005

Your Local Agent

QLD 6-8 Imboon St, Deception Bay 4508
Ph 1300 368 085 qld@rehabandmobility.com.au
www.rehabandmobility.com.au