



Bathroom

Aluminium Shower Chair RM402

The RM lightweight aluminium shower chair with durable “S” bend side frame offering maximum comfort and space when seated. With large rubber suction ferrule to prevent slipping and provide maximum safety. Moulded plastic seat with drainage holes.



- * Aluminium frame
- * Plastic seat with drainage holes
- * “S” bend side frame provides maximum space
- * Large rubber suction ferrule to prevent slipping
- * Weight capacity 160kg

Specifications

Overall height.....	79-94	Seat height	44-59
Overall depth	64	Armrest height	21
Overall width.....	60	Unit weight.....	3.3
Seat depth.....	41	Weight capacity	160
Seat width	49		

Measurements in cm, weight in kg

Product Code: RM402

Your Local Agent

QLD 6-8 Imboon St, Deception Bay 4508
Ph 1300 368 085 qld@rehabandmobility.com.au
www.rehabandmobility.com.au