

Tubular Slide Sheets

The Aidapt range of slide sheets are manufactured from a high quality silicone coated nylon. The low friction slide sheets aid and assist carers in the safe moving and transferring of patients. These slide sheets reduce the risk to the carer from strains and injury. The heat sealed edges are ideal where there is a high risk of cross contamination between patients. The tubular slide sheets come in a variety of sizes.



Product Code: VM500



Product Code: VM500B



Product Code: VM500Y



Product Code: VM500G

Specifications

VM500 Red 600X400X3mm

VM500B Blue 720x700x3mm

VM500Y Yellow 1450x710x3mm

VM500G Green 1220x1000x3mm

- * Tubular formed for easy positioning of a patient
- * Diminishes force being applied to patient
- * Low friction
- * Made of nylon coated in silicone

Glide Sheet with Handles



- * Easy transferring of patient
- * Reduces strain on patient and carer
- * Made from nylon coated in silicone
- * Comes with handles

Specifications

Dimensions 1900x1000mm

Product Code: VM500R

Your Local Agent

QLD 6-8 Imboon St, Deception Bay 4508
Ph 1300 368 085 wholesale@rehabandmobility.com.au
NSW Unit 2/144 Lake Entrance Road, Oak Flats 2529
Ph 1300 002 805 nsw@rehabandmobility.com.au
www.rehabandmobility.com.au