



Fact Sheet

Bathroom



Solo Wall-mounted Seats

Solo wall-mounted stools and chairs provide a compact solution, offering a secure showering seat that can be folded away against the wall when not in use, ideal for areas where space is limited.



Product Code: VB544

- * All alloy lightweight frame, robust construction
- * Anti-bacterial finish to the ergonomic contoured padded PU seat
- * Fine adjustment leg for extra stability

Specifications

Seat dimensions	36x28
Seat height	47-54
Weight capacity	159



Product Code: VB634

- * Lightweight aluminium
- * Simple to install and maintain
- * Clip on, clip off seat with drainage for easy cleaning

Specifications

Seat dimensions	49x38
Seat height -stool	46-54
Seat height -chair	48-61
Width between arms.....	46
Depth from wall	49
Weight capacity.....	190

Measurements in cm, weight in kg



Product Code: V535

Your Local Agent

QLD 6-8 Imboon St, Deception Bay 4508
Ph 1300 368 085 qld@rehabandmobility.com.au
www.rehabandmobility.com.au