



Bathroom



Multi Purpose Stool with Rotating Seat

This versatile, multi-purpose stool can be used in a whole host of different settings around the home. With a comfortable rotating seat top making it ideal for those who need assistance when turning from side to side when seated.



- * Durable, lightweight rust proof aluminium frame
- * Suitable for use in bath / bathroom
- * Rotating seat top
- * Height Adjustable, 8 settings
- * Suitable for children or adults
- * Non-slip rubber tips
- * Retail packaging
- * Weight capacity 135kg



Specifications

Overall height 35-52
 Overall depth 33
 Overall width 33
 Seat depth 33
 Seat width 33

Seat height 35-52
 Unit weight 2.5
 Weight capacity 135

Measurements in cm, weight in kg

Product Code: VB511R

Your Local Agent

QLD 6-8 Imboon St, Deception Bay 4508
 Ph 1300 368 085 qld@rehabandmobility.com.au
 www.rehabandmobility.com.au