



Fact Sheet

Bathroom



Toilet Surround

Height adjustable Toilet Surrounds assist in getting on and off the toilet without having to install wall rails. Durable aluminium powder coat finish.

- * Height adjustable
- * Width adjustable (VR202 & VR205SP)
- * Double transform bar for added strength
- * Coated with a tough white polymer anti-scratch finish
- * Moulded support arms provide improved grip
- * The toilet frame dismantles easily for storage/transportation
- * Will fit most toilet suites
- * Rubber tips to reduce slips



Product Code: VR202



Product Code: VR205SP

- * Detachable frame
- * Available in 54cm & 60cm width



Product Code: RM418-54 (54cm)
RM418-60 (60cm)

Specifications VR202 & VR205

Overall height 66-81
Overall depth 45
Overall width 58
VR202 int width 49-63
VR205 int width.....47.5-52.5

Unit weight 4.3
Weight capacity 190

Measurements in cm, weight in kg

Specifications RM418

Overall height 65-77
Overall depth 52
Overall width 58/62
Unit weight 2.6
Weight capacity 130

Your Local Agent

QLD 6-8 Imboon St, Deception Bay 4508
Ph 1300 368 085 qld@rehabandmobility.com.au
www.rehabandmobility.com.au