

Fact Sheet



Seating Bariatric



Torbay Perching Stool

The Torbay Bariatric Perching Stool has been developed for use in the kitchen and around the home, suitable for those who find it difficult to stand for longer periods. With larger footprint for stability, adjustable height, wide seat frame and reduced incline of seat for comfort and safety. Higher arm supports to aid getting on and off the stool.

- * Ideal for perching when working in the kitchen or at a bench
- * Adjustable seat height 7 settings
- * Easy clean
- * Weight Capacity 254kg



Product Code: VG868W



Product Code: VG869W



Product Code: VG871W

Specifications—VG871W

Overall height..... 96-101
Overall depth51
Overall width.....69
Seat depth29
Seat width48

Seat height 54-71
Armrest height23
Unit weight.....9.1
Weight capacity.....254

Measurements in cm, weight in kg

Your Local Agent

QLD 6-8 Imboon St, Deception Bay 4508
Ph 1300 368 085 qld@rehabandmobility.com.au
www.rehabandmobility.com.au