

# Fact Sheet



## Bathroom

### Tri Corner Stool with Padded Seat

The Tri Corner Shower Stool is an excellent option when space is at a minimum. Lightweight, sturdy with rubber PVC stoppers for anti-slip stability. Contoured padded seat makes this chair a comfortable stool for use in the bathroom.



- \* Ideal where space is limited
- \* Lightweight aluminium (3.16kg)
- \* Height adjustable
- \* Drain holes under padded seat
- \* Non-slip rubber tips
- \* Weight capacity 136kg

#### Specifications

Overall height ..... 49-65  
 Overall depth ..... 45  
 Overall width ..... 42  
 Seat depth ..... 45  
 Seat width ..... 42

Seat height ..... 49-65  
 Unit weight ..... 3.1  
 Weight capacity ..... 136

*Measurements in cm, weight in kg*

**Product Code: RM396PU**

**Your Local Agent**

**QLD** 6-8 Imboon St, Deception Bay 4508  
 Ph 1300 368 085 [qld@rehabandmobility.com.au](mailto:qld@rehabandmobility.com.au)  
[www.rehabandmobility.com.au](http://www.rehabandmobility.com.au)