

Fact Sheet



Walking Aids



Tri Walker

The aluminium Tri-Walker from Aidapt is sturdy yet easy to manoeuvre, featuring lockdown brakes. Suitable for both indoor and outdoor use, supplied with bag and basket with tray easily folds away for transport and storage.

- * Ideal for indoor and outdoor use
- * Adjustable height handles for a comfortable fit
- * Large 8" castor
- * Easy folding for storage and transport
- * Simple to use lock brakes
- * Supplied with bag, basket and tray
- * Available; silver, blue & red



Specifications

Overall height83-99
 Overall depth 60
 Overall width 67
 Folded width 28

Unit weight 6.8
 Weight capacity 133

Measurements in cm, weight in kg

Product Code: VP172AA Red

Product Code: VP172AB Blue

Product Code: VP172AC Silver

Your Local Agent

QLD 6-8 Imboon St, Deception Bay 4508
Ph 1300 368 085 qld@rehabandmobility.com.au
www.rehabandmobility.com.au