



Fact Sheet

Wheelchair



Urban Wheelchair

Quality Foldable Transit Wheelchair makes life a bit simpler, for both the occupant and the attendant. A smart, sporty and classy looking compact unit, it is comfortable and user friendly.

- * Ultra light frame
- * Folding backrest
- * Attendant Brake
- * Swing up leg rests
- * Rear wheel locking brake
- * Weight capacity 120kg



Specifications AP/SP

Overall height	86	Seat height	46
Overall depth.....	100/112	Backrest height.....	41
Overall width.....	58/64	Armrest height	23
Seat depth	40	Unit weight.....	10.5/11.5
Seat width	46	Weight capacity.....	120

Measurements in cm, weight

Product Code: WCT16 A/P
Product Code: WCT20 S/P

Your Local Agent

QLD 6-8 Imboon St, Deception Bay 4508
Ph 1300 368 085 qld@rehabandmobility.com.au
www.rehabandmobility.com.au