



Fact Sheet

Pressure Care



Viscotech Contour Plus

Viscotech® Plus Contour is a specially contoured composite visco-elastic and foam cushion providing excellent pressure relieving qualities and superb comfort. The welded cover offers maximum infection prevention with two way stretch and non-slip base. Viscotech® Plus Contour offers positioning, is suitable for the active user and gives supreme envelopment to the client.



- * High Risk
- * Contoured composite visco-elastic and foam
- * Waterproof, vapour permeable polyurethane cover
- * Welded seams to prevent fluid entry
- * Wide flat covered zip
- * Made in Britain
- * Weight Capacity 108kg

Specifications

Product Code: VCC1816W 46x41x9cm
VCC2018W 51x46x9cm

Your Local Agent

QLD 6-8 Imboon St, Deception Bay 4508
Ph 1300 368 085 wholesale@rehabandmobility.com.au
NSW Unit 2/144 Lake Entrance Road, Oak Flats 2529
Ph 1300 002 805 nsw@rehabandmobility.com.au
www.rehabandmobility.com.au