

Walking Aids



Walking Dual Riser Frame

A Dual Rise Lightweight Aluminium Folding Walking Frame with a wide base for maximum stability, ideal for helping users to rise safely from a seated position. The lightweight aluminium walking frame features soft grip handles, one button folding function and height adjustable legs. Can be supplied with wheels.



Product Code: VP179A

- * Lightweight Aluminium frame
- * Dual handle for easy rising
- * Riser handle distance from top handle 290mm
- * Easy push button fold
- * Non-slip rubber stoppers
- * Can be fitted with wheels or skis
- * Max weight user: 133kg

Walking Frame Accessories

VP179T Apron Bag

VP179N Net Bag

(details available www.rehabandmobility.com.au)

Specifications

Overall Height	Overall Length	Overall Width	Folded Width	Unit Weight	Weight Capacity
82-99	52	61	12	2.6	133

Measurements in cm, Weight in kg

Your Local Agent

QLD 6-8 Imboon St, Deception Bay 4508
Ph 1300 368 085 wholesale@rehabandmobility.com.au
NSW Unit 2/144 Lake Entrance Road, Oak Flats 2529
Ph 1300 002 805 nsw@rehabandmobility.com.au
www.rehabandmobility.com.au