



# Fact Sheet

## Walking Aids



### Walking Dual Riser Frame

A Dual Rise Lightweight Aluminium Folding Walking Frame with a wide base for maximum stability, ideal for helping users to rise safely from a seated position. The lightweight aluminium walking frame features soft grip handles, one button folding function and height adjustable legs. Can be supplied with wheels.



- \* Lightweight Aluminium frame
- \* Dual handle for easy rising
- \* Riser handle distance from top handle 290mm
- \* Easy push button fold
- \* Non-slip rubber stoppers
- \* Can be fitted with wheels or skis
- \* Max weight user: 133kg

#### Walking Frame Accessories

VP179T Apron Bag

VP179N Net Bag

(details available [www.rehabandmobilityshop.com.au](http://www.rehabandmobilityshop.com.au))

#### Specifications

Overall height..... 82-99  
 Overall depth .....52  
 Overall width.....61  
 Folded width .....12

Unit weight.....2.6  
 Weight capacity.....133

*Measurements in cm, weight in kg*

**Product Code: VP179A**

#### Your Local Agent

**QLD** 6-8 Imboon St, Deception Bay 4508  
**Ph 1300 368 085** [qld@rehabandmobility.com.au](mailto:qld@rehabandmobility.com.au)  
[www.rehabandmobility.com.au](http://www.rehabandmobility.com.au)